 **esperanza**[®]
creating space for your dreams

BABY SCALE

BEBE

MANUAL
EBS019

BABY SCALE

IMPORTANT SAFETY INSTRUCTIONS

- Carefully unpack the device and take care not to throw away any part of the packaging material before you have found all parts of the product.
- Under no circumstances should you repair or alter the product by yourself. Entrust all repairs and adjustments to a specialized company/service. Failure to do so may void the warranty.
- Store and use the product away from the flammable or volatile materials and substances.
- Do not expose the product to extreme heat, direct sunlight or excessive humidity and do NOT place in an overly dusty environment.
- Do not use this product for other purposes than those intended.
- Protect the device from impact with hard objects, moisture, dust, chemicals, toiletries, liquid cosmetics, great temperature fluctuation and closeness to sources of heat (open fires, radiators).
- If you will not be using this product for an extended period of time, remove the batteries. Otherwise they may leak inside the product and cause irreparable damage.
- Keep far away from strong electromagnetic fields (e.g. cell phones); This may impair the accuracy of the scale.

CHILDREN AND FRAIL INDIVIDUALS

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach. WARNING! Do not allow small children to play with the foil as there is a danger of suffocation!
- This device is not intended to be used by individuals (including children) who have restricted physical, sensory or mental abilities and/or insufficient knowledge and/or experience, unless they are supervised by an individual who is responsible for their safety or have received instructions on how to use the device.
- Children should be supervised at all times in order to ensure that they do not play with the device.

SPECIFICATIONS

- Precise weighing accuracy
- Capacity: 0,03- 20 kg
- Scale division: 0,01kg
- Tare function
- Automatic hold function
- Big LCD display with green back light: 72*27 mm
- Unit: KG/LB/OZ
- 150 cm soft measuring tape with automatic roll back function
- Power supply: 2* AAA Battery (not included)
- Low battery and overload indication
- Automatic and manual power off
- work temperature: 0°C - 40°C
- Product size: 54,5 *32*3 cm

OPERATION

Weighing and hold function

- Place the scale on a level and solid surface. On the back of scale is the battery compartment. Insert batteries into compartment following the polarity direction which is marked inside the compartment. Put the cover back.

- Press the button "ON/OFF/TARE". When the scale is placed on a not level or solid surface, the LCD will show "UnST". Switch off the scale and place it on a level and solid surface. Press the button "ON/OFF/TARE" again, LCD shows "0.00", the scale is ready to take measurements.
- Press the button "UNIT/ HOLD" shortly, you can change the unit between "kg", "lb", "oz". Select your desired unit, LCD shows "0.00" together with the unit you selected. Hold and press the button "UNIT/ HOLD" for 4 seconds, LCD display " L_ON" which means the function of automatic holding the measurement result is switched on. Hold and press the button "UNIT/ HOLD" for 4 seconds again, LCD display " L_OF" which means the function of automatic holding the measurement result is switched off.
- Place baby on the platform of scale, the weight is now displayed. When baby moves a lot, the weight is also changing fast. When baby is more calm and weight changes getting slower and slower, in case the automatic holding function is switched on, the weighing result will flash 2 times and then it is blocked on the display with the " hold" displayed on the left top of LCD.

Tare Function

If you do not want to place the baby directly on the weighing surface, but prefer to use a towel and/or any other items which you want to deduct during weighing, you can use tare function.

- Press button "ON/OFF/TARE" to switch on the scale. After LCD shows "0.00", place an towel on the weighing surface, in case of the automatic holding function is on, the weight of towel will flash 2 times. Press button "ON/OFF/TARE", LCD shows: "TArE" and then "0.00", Place baby on the scale. When the weight getting more and more stable, the weighing result will flash 2 times. Now it shows the weight of baby. Press button "ON/OFF/TARE" again, LCD shows: "TArE" and then "0.00". When you remove the towel or blanket, the "- weight" of the towel will appear on the screen, because the weight data is negative now.
- Press button "ON/OFF/TARE" again, "- weight" will disappear and then "0.00" shows up on screen. You can continue weighing or continue clear other items.

MEASURING THE BABY HEIGHT

On the right side of the scale, there is a soft measuring tape with a total length 150cm. Easily pull it out, you can measure your baby height.

When you finish measuring, just press the button on the back of scale, the soft measuring tape will roll back.

Switching off

Press the key of "ON/OFF/TARE" for 3 seconds, the scale will be power off. Or the scale will power off automatically after 30 seconds when there is no action during this time.

ERROR DISPLAY INTRODUCTION

- When the scale is over loaded, the "O-Ld" will appear in the display.
- When "Lo" appears in the display, the batteries need to be replaced.
- If "- weight" show on the screen, that means weight result is negative. Press "ON/OFF/TARE" button to re-activate the scale.

MAINTAINANCE AND STORAGE

- Cleaning: The scale should be cleaned with a damp cloth and a little washing up liquid if required.
- Never immerse the scale in water or rinse it under running water.
- Storage: Do not place any objects on the scale when not in use.

WARRANTY

1. This warranty is issued by Esperanza sp.j. Poterek with address:
ul. Poznańska 129/133, 05-850 Ożarów Mazowiecki, Phone: 022-7213571 w.10 – called hereinafter Esperanza.
2. This warranty confers upon you specific rights and remedies against Esperanza. Please read this document carefully, by making a claim under this warranty, you will be deemed to have understood and accept its terms.
3. Esperanza warrants that any Esperanza product together with attached accessories is free from defects in materials and workmanship under conditions of its normal use.
4. If you make a claim under this warranty, Esperanza will, at its option, repair or replace any defective product into new one by using new or use factory-refurbished parts. In case of using factory-refurbished parts they will be equal to new products.
5. Defective product should be returned together with receipt, commercial invoice or any other proof of purchase to the place of purchase. It should include all accessories attached in original packing. Esperanza has right to reject warranty if any accessories will be missed.
6. Esperanza is committed to respond to the complaint within 14 calendar days of its notification in the point of sale. Esperanza will repair the defective product in the shortest possible time. If will be need to import from abroad of any parts which are necessary to repair this time may be extended.
7. In case of replacing defective product onto new one, Esperanza has right to replace it with newer model with the same technical parameters if the old one model was replaced by new one during warranty period.
8. If you make claim under this warranty, you are deemed to have agreed that repair and (upon availability) replacement, as applicable, under the warranty services described herein, is your sole and exclusive remedy vis-a-vis Esperanza.

EXCLUSIONS:

This warranty does not extend to:

1. Any product not manufactured by Esperanza
2. Any product that has been damaged or rendered defective as a result of:
 - use it in other way than normal intended use
 - modification or internal interferes or changes of product
 - service or repair by anyone other than Esperanza service
 - improper store, transportation or packing
 - improper installation of product or software
 - mechanical, physical, chemical etc. damages
3. In no event shall Esperanza be liable for any damages, loss of business, profits, data or loss based on use Esperanza products with other devices.
4. Products with removed warranty and other labels and products spent or have hard use marks.

All products are covered by Esperanza's 2-year manufacturer's warranty while product is purchased by the consumer (the consumer is a person who buys product for the purpose of not professional or commercial activity). In the case of purchase for the company (for purposes relating to the activities of professional or commercial) warranty is 1 year from date of purchase. Esperanza products have all the necessary safety certificates required in the EU.

According to Directive 2012/19/EU, disposal of waste electrical and electronic equipment together with household waste is prohibited. Their components need to be recycled or disposed separately. Otherwise, they can cause the release of hazardous and noxious substances which pose a threat to the environment and human health. The law requires each consumer to return unnecessary or waste electrical and electronic equipment to the producers, sales agents or waste treatment facilities. Detailed conditions are governed by national legislation. The above symbol placed on the product, in the user manual or on the packaging refers to this requirement. Thanks to proper segregation and compliance with the provisions relating to waste disposal each user contributes significantly to environmental protection.

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